LEVELON-VR

- Levelling Agent For Fluorescent Dyestuff

*** SPECIFICATION**

* APPEARANCE

* IONIC NATURE

* pH (1% aq. Sol'n)

* SP.GR(at 25°C)

* SOLUBILITY

* STORAGE STABILITY

* COMPATABILITY

Yellowish Liquid

Anionic

 8.0 ± 1.0

 1.00 ± 0.05

Easily soluble in warm and cold water.

Min. 1 year under cool and dry conditions when stored in original sealed, packaging.

Compatible with anionic and nonionic materials.

*** FEATURES**

- 1) In dyeing nylon and wool with fluorescent dyestuff, **LEVELON-VR** gives good migration and retardation effects of dyestuff.
- 2) LEVELON-VR provides good dyeing affinity in dyeing with dyestuffs with poor dyeing affinity such as Turquoise-Blue and Fluorescent dyestuffs as well as acid dyestuff.
- 3) Clearer fluorescent color can be obtained when **LEVELON-VR** is used as a levelling agent in dyeing with fluorescent dyestuff.
- 4) **LEVELON-VR** provides excellent levelling effect by improving penetration of dyestuff.
- 5) **LEVELON-VR** cuts down the cost because it can save approx. 20% dyestuff by improving dyeing affinity of dyestuff with poor dyeing affinity like Turquoise-Blue dyestuff.
- 6) **LEVELON-VR** prevents tarring of dyestuff in dyeing and has excellent migration effect.
- 7) It is environmental-friendly product as APEO-free and NPEO-free. Suitable for Nylon and Wool fibers

* APPLICATION FIBERS

- **HOW TO USE**
- RECOMMENDEDRECIPE

Recommended use levels are 1-2% (o.w.f) according to kinds of dyed materials, kinds of dyestuffs, concentration of dyestuffs and dyeing methods.


Acid Dyestuff : X % (O.W.F.)
 LEVELON-VR : 1 ~ 2% (O.W.F.)

• Acetic acid : Adjust to the pH $3.5 \sim 4.0$

• Bath Ratio : 1:30 ~ 1:50

• Time and Temp. : $60 \text{ min at } 95 \sim 100^{\circ}\text{C}$

DYEING DIAGRAM



CORRECTION METHOD OF DYEING SPOTS


- Modification of the faulty dyeing of wool and knitted nylon fabrics is as follows: use 2~3% LEVELON-VR (O.W.F.), and treat it for 30~45 minutes at 85~95°C, pH 3~4.
- In case it is used with fluorescent dyestuff that has good levelling effect,
 LEVELON-VR provides same effect even with 50% of the required amounts.
- If the faulty dyeing occurs although goods are dyed as above method, use
 LEVELON-VR according to the correction process of the faulty dyeing as below.

• **LEVELON-VR** : 2 ~ 3% (O.W.F.)

• Acetic acid : Adjust to the pH $3.5 \sim 4.0$

• Bath Ratio : 1: 30 ~ 1: 50

• Time and Temp. : 30 min at $95 \sim 100^{\circ}$ C



PACKAGE

120 Kgs NET IN DRUMS.